dentle | zip

USER MANUAL INFORMATION

BRUKSANVISNING

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KÄYTTÖOPAS

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CONTENTS

- Electric Toothbrush (1)
- Brush Head (2), On/Off Button (3), Battery Cover (4), Case (5)

TECHNICAL DATA

- Battery: 1 x 1.5V AAA size alkaline battery
- Dimensions: width 22.2mm, length 221.5mm (with cover)
- Weight: 40.5g (with cover and battery)
- Main unit: IPX4

SPECIFICATIONS

- Vibrates 23,000 times per minute
- A new alkaline battery lasts for up to three months when brushing twice a day for two minutes each time

INSTRUCTIONS

Assembly

1. The toothbrush comes in two parts that are easy to assemble.
2. Simply insert the toothbrush into the cover and it becomes a normal-sized toothbrush.

Brushing your teeth

1. Begin by dipping the brush head into water. Then squeeze out a small amount of toothpaste onto the brush head.
2. Avoid splashing by waiting until the toothbrush is in contact with your teeth before turning it on.
3. Brush your teeth for about two minutes. Move the brush head from tooth to tooth. Begin by brushing the outside surfaces, then the inside ones and, finally, the chewing surfaces.
4. Rinse the brush head when you have finished and let it dry.

Changing the brush head

1. Remove the old brush head by turning the brush head in an anti-clockwise direction and then pulling it straight up.
2. Attach the new brush head by pressing it down and then turning the brush head in a clockwise direction.

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