# QUICK GUIDE

on how to operate your Dente One Plus toothbrush

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## SAFETY RULES

Read the Quick Guide before using the toothbrush. • Never handle the charger adapter with wet hands. • Never put the charging plate into water or any other liquid. • Never force the charger adapter into an electrical socket. • Never allow a child to use the Dentle One without adult supervision. • Never use the Dentle One while taking a bath or shower. • Make sure that the cord does not come into contact with hot surfaces. • Never use the Dentle One if the charger adapter is damaged. • Dispose of the unit if the charger adapter is damaged. The charger adapter cannot be replaced. • Use the Dentle One only for its intended purpose, as described in the Quick Guide. This electric toothbrush is designed only for cleaning teeth, and not for any other use. • Use only accessories that are recommended by the manufacturer. • Replace the brush head if the bristles are worn out to avoid having the bristles fall off or loosen when the electric toothbrush is being operated. • Never use the charger if it has accidentally been dropped in water. • Never bend or force the power plug into the wrong type of socket or power strip. • Never disassemble the Dentle One yourself. This should only be done by specialised personnel. • If you experience any irritation or adverse reactions while using the Dentle One, please discontinue its use and consult your dentist. • Children should be supervised to ensure that they do not play with the appliance. • This appliance can be used by children aged from 8 years and above and persons with reduced physical, sensory or mental capabilities or lack of experience and knowledge if they have been given supervision or instruction concerning use of the appliance in a safe way and understand the hazards involved. Children must not play with the appliance. Cleaning and user maintenance must not be performed by children without supervision. • For hygiene reasons the charging cup should only be used for its intended purpose of charging the electric toothbrush. It should neither be used for drinking nor be filled with water or other liquids. • Only use the toothbrush with the provided power supply and charging station.

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**dentle** | one plus is an electric toothbrush designed in Sweden. A product's appearance is an important part of the total experience, so every part of the Dentle One Plus combines elegant simplicity with superior performance and dependability. The electric toothbrush, the plug and the charging station cup, have all been designed to be easy to use and attractive to look at, thus contributing to healthier teeth and gums, while fitting nicely into your home. This Quick Guide contains simple instructions about how to use the Dentle One Plus. At [dentle.com](http://dentle.com), you will find additional advice and suggestions for how to best care for your teeth, produced in cooperation with Swedish dentists. We are pleased that you chose to purchase a Dentle One Plus, and hope you'll be very satisfied with your new electric toothbrush.

## CONTENTS

Electric toothbrush **A**, brush head **B**, on/off button **C**, charge indicator **D**, charging cup (supplied in three parts) **E**, charging adapter with cord and docking plug **F**.

### Technical data

Battery: 1 x ICR14500 800 mAh Li-ion battery.  
Dimensions: width 30 mm, length 222 mm. Weight: 230 g (electric toothbrush and charger).  
Adaptor: I/P: 100–240 V AC; 50/60 Hz; 0.2 A; Class II; IPX7 O/P: 5 V; 1A. Charging base: 5 V, 1A; Class III; IPX7 except for the jack socket. Toothbrush: 3.7 V; Class III; IPX7 I/P: 1.8 W.

## SPECIFICATIONS

- ⌚ Oscillates clockwise and counter-clockwise 10,000 times per minute.
- ⌚ Three speed settings: high (10 000 oscillations), medium (8 500 oscillations) and low (7 000 oscillations). Press and hold the power button for 1 second to toggle mode.

- ⌚ Timer function that automatically switches the brush off after two minutes, and indicates after 30, 60 and 90 seconds that it is time to brush a different area.
- ⌚ Charge lasts up to 28 days when brushing twice a day for two minutes each time.
- ⌚ Charge initially for 14 hours.
- ⌚ The charging time for the toothbrush depends on how often you charge it. For the best charging time, charge your toothbrush every time you use it.
- ⌚ In order to get the best possible results, and thus healthier teeth and gums, remember to change the brush head every month.

## INSTRUCTIONS

### Assembly

- 1 Attach the brush head by pressing it down with the bristles facing the same way as the button on the handle.
- 2 The charging cup comes in three parts that are easy to assemble. The smaller ring part clicks onto the charging plate. Attach the larger ring to the smaller one to complete assembling the cup.
- 3 Connect the cord of the charging adapter to the cup and then connect the adapter to an electrical outlet.
- 4 Place your electric toothbrush in the charging cup to begin charging your Dentle One. The charge indicator flashes red when the battery needs charging, pulsates white as soon as the battery is at 20% charge and shows a steady white light when fully charged.
- 5 Your Dentle One electric toothbrush should be charged for at least 14 hours before first use.
- 6 Disconnect the charger as soon as the electric toothbrush is fully charged.

### Brushing your teeth

- 1 Begin by dipping the brush head into water. Then squeeze a little toothpaste onto the toothbrush head.
- 2 Avoid splashing by waiting to turn the toothbrush on until it is in contact with your teeth.
- 3 You can choose between three different speed modes. Rotate through the modes by pressing the on/off button, each time keeping it down for 1 second. Choose between: high (10,000 oscillations), medium (8,500 oscillations) or low (7,000 oscillations).
- 4 Brush your teeth for approx. 2 minutes (the timer function automatically switches the brush off after 2 minutes). Move the brush head from tooth to tooth. Start with the outside of the teeth, then move on to the inside and finish with the chewing surfaces. The timer function indicates after 30, 60 and 90 seconds that it is time to start brushing the next quadrant of the mouth.
- 5 Rinse off the brush head when you have finished, and allow it to dry.

### Changing the brush head

- 1 Remove the old brush head by pulling it straight up.
- 2 Attach the new one by pressing it down with the bristles facing towards you.
- 3 It is important to change the brush head once a month for best results.

### Dismantling the charging cup

It is possible to dismantle the charging cup when necessary:

- 1 Lift off the larger ring from the smaller ring.
- 2 Then lift the smaller ring from the charging plate, and that's it.

## RECYCLING

The toothbrush contains a rechargeable lithium-ion battery. The battery is recyclable and can be left at a recycling station.

### To remove the battery for disposal or recycling:

- 1 The appliance must be disconnected from the supply mains when removing the battery.
- 2 Insert a coin into one set of slots on the battery end cap and twist counterclockwise to release the battery cap.
- 3 Firmly tap the bottom of the handle on a hard surface in order to slide out battery.
- 4 Now the battery is exposed so you can separate the wire connectors.
- 5 Disconnect the individual cables one after another and remove the rechargeable batteries.
- 6 The battery must be removed from the appliance before it is scrapped.
- 7 Safely and properly dispose of the battery.

### Correct Disposal of this product

 This marking indicates that this product should not be disposed with other household wastes throughout the EU. To prevent possible harm to the environment or human health from uncontrolled waste disposal, recycle it responsibly to promote the sustainable reuse of material resources. To return your used device, please use the return and collection systems or contact the retailer where the product was purchased. They can take this product for environmental safe recycling.

The Quick Guide is also available at Dentle's homepage [www.dentle.com](http://www.dentle.com).

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Super Living Products Ltd. Room 408B, Lippo Sun Plaza, 28 Canton Road, Kowloon, Hong Kong Model No. 7806