QUICK GUIDE

on how to operate your Dentle One toothbrush

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SAFETY RULES

- Read the Quick Guide before using the toothbrush.
- Never handle the charger adapter with wet hands.
- Never put the charging plate into water or any other liquid.
- Never force the charger adapter into an electrical socket.
- Never allow a child to use the Dentle One without the supervision of an adult.
- Never use the Dentle One when taking a bath or shower.
- Make sure that the cord does not come in contact with hot surfaces.
- Never use the Dentle One if the charger adapter is damaged.
- Dispose of the unit if the charger adapter is damaged. The charger adapter cannot be replaced.
- Use the Dentle One only for its intended purpose, as described in the Quick Guide. This electric toothbrush is designed only for cleaning teeth, and not for any other use.
- Use only accessories that are recommended by the manufacturer.
- Replace the brush head if the bristles are worn out to avoid having the bristles fall off or loosen when the electric toothbrush is being operated.
- Never use the charger if it has accidentally been dropped in water.
- Never bend & force the power plug into the wrong type of socket or power strip.
- Never disassemble the Dentle One by yourself. This should only be done by specialized personnel.
- If you experience any irritating or adverse reactions while using the Dentle One, please discontinue its use and consult your dentist.
- Children should be supervised to ensure that they do not play with the appliance.
- This appliance is not intended for use by persons (including children) with reduced physical, sensory or mental capabilities, or lack of experience and knowledge, unless they have been given supervision or instruction concerning use of the appliance by a person responsible for their safety. Cleaning and user maintenance shall not be made by children without supervision.
- For hygienic reasons the charging cup should only be used for its intended purpose of charging the electric toothbrush. It should neither be used for drinking out of nor be filled with water or other liquids.

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dentle |one is an electric toothbrush designed in Sweden. A product's appearance is an important part of the total experience, so every part of the Dentle One combines elegant simplicity with superior performance and dependability. The electric toothbrush, the plug and the charging station cup, have all been designed to be easy to use and attractive to look at, thus contributing to healthier teeth and gums, while fitting nicely into your home. This Quick Guide contains simple instructions about how to use the Dentle One. At dentle.com, you will find additional advice and suggestions for how to best care for your teeth, produced in cooperation with Swedish dentists. We are pleased that you chose to purchase a Dentle One, and hope you'll be very satisfied with your new electric toothbrush.

CONTENTS

The Dentle One introduction package contains a rechargeable electric toothbrush, a brush head, a charging cup (delivered in three parts), a charger adapter with a cord and docking plug to the charging cup.

Contents

Electric toothbrush ①, brush head ②, on/off button ③, charge indicator ④, charging cup (delivered in three parts) ⑤, charging adapter ⑥ with cord and docking plug ⑦.

Technical data

Batteries: 2 x 3/5 AA 600mAh NiMH batteries.
Dimensions: width 30 mm, length 222 mm. Weight: 230 g (electric toothbrush and charger). Adaptor: I/P: 100-240VAC; 50/60Hz; 1,8W; Class II; IPX7 O/P: 26VAC; 120mA. Charging base: 26VAC; Class III; IPX4. Toothbrush: 2.4V; Class III; IPX7 Input power 1.6W.

SPECIFICATIONS

- ④ Oscillates clockwise and counter-clockwise 8,000 times per minute.
- ④ Three speed settings: high (8,000 oscillations), low (7,000 oscillations) and massage (6,000–8,000 oscillations).
- ④ Timer function that automatically switches the brush off after two minutes, and indicates after 30, 60 and 90 seconds that it is time to brush a different area.
- ④ Charge lasts up to 14 days when brushing twice a day for two minutes each time.
- ④ Charge initially for 24 hours.
- ④ The charging time for the toothbrush depends on how often you charge it. For the best charging time, charge your toothbrush every time you use it.
- ④ In order to get the best possible results, and thus healthier teeth and gums, remember to change the brush head every month.

INSTRUCTIONS

Assembly

- ① The charging cup comes as three parts that are easy to assemble. The smaller ring part clicks onto the charging plate. Just attach the larger ring to the smaller one, and you've finished.
- ② Connect the charging adapter to an electric outlet. Then place your electric toothbrush in the charging cup, and your Dentle One will start charging.
- ③ Your Dentle One electric toothbrush should be charged for at least 24 hours before using it for the first time.

- ④ The charging indicator (4) lights up red when the battery needs charging, flashes green while the battery is charging, and lights up constantly green when the battery is fully charged.
- ⑤ Disconnect the charger as soon as the electric toothbrush is fully charged.

Brushing your teeth

- ① Begin by dipping the brush head into water. Then squeeze out a bit of toothpaste onto the toothbrush head.
- ② Avoid splashing by waiting to turn on the toothbrush until it is in contact with your teeth.
- ③ Brush your teeth for approx. 2 minutes (the timer function automatically switches the brush off after 2 minutes). Move the brush head from tooth to tooth. Start with the outside of the teeth, then move on to the inside and finish up with the chewing surfaces. The timer function indicates after 30, 60 and 90 seconds that it is time to start brushing a different area of the mouth.
- ④ Rinse off the brush head when you are finished, and let it dry.

Changing the brush head

- ① Remove the old brush head by pulling it straight up.
- ② Attach the new one by pressing it down with the bristles facing towards you.
- ③ It is important to change the brush head once a month for the best results.

Dismantling the charging cup

It is possible to dismantle the charging cup when necessary.

- ① Lift off the larger ring from the smaller ring.
- ② Then lift the smaller ring from the charging plate, and you're done.

RECYCLING

This toothbrush contains two rechargeable nickel-metal hydride batteries (NiMH). These batteries are recyclable, and can be left at a recycling station for recycling.

To remove the battery for disposal or recycling:

- ① The appliance must be disconnected from the supply mains when removing the battery.
- ② Insert a coin into one set of slots on the battery end cap and twist counterclockwise to release the battery cap.
- ③ Firmly tap the bottom of the handle on a hard surface in order to slide out battery.
- ④ Now the battery is exposed so you can separate the wire connectors.
- ⑤ Disconnect the individual cables one after another and remove the rechargeable batteries.
- ⑥ The battery must be removed from the appliance before it is scrapped.
- ⑦ Safely and properly dispose of the battery.

Correct Disposal of this product

 This marking indicates that this product should not be disposed with other household wastes throughout the EU. To prevent possible harm to the environment or human health from uncontrolled waste disposal, recycle it responsibly to promote the sustainable reuse of material resources. To return your used device, please use the return and collection systems or contact the retailer where the product was purchased. They can take this product for environmental safe recycling.

The Quick Guide is also available at Dentle's homepage www.dentle.com.

The charging cup comes as three parts that are easy to assemble.

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